



## Product Specification

**S-Monovette® Lithium heparin LH, 9 ml, cap orange, (LxØ): 92 x 16 mm, with plastic label**

**Page 1**



### Product description

Order number	02.1065
Product description	S-Monovette® Lithium heparin LH, LH, preparation: Lithium heparin, 9 ml, membrane screw cap, cap orange, colour code EU, (LxØ) without cap: 92 x 16 mm, with plastic label, label/print: transparent/orange, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Lithium heparin
Preparation concentration	16 IU/ml blood
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	transparent/orange
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® Lithium heparin LH, 9 ml, cap orange, (LxØ): 92 x 16 mm, with plastic label**

Page 2

Liquid-tight Connection yes  
membrane screw cap

### Size

Sample volume 9 ml  
Outer diameter of cap 19 mm  
Diameter 16 mm  
Length including cap 108 mm  
Length excluding cap 92 mm  
Length including cap and plunger 129 mm

### Material & colours

Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap orange  
Membrane material rubber  
Tube material Polypropylene (PP)  
Piston material High Density Polyethylene (HD-PE)  
Piston plunger material Polystyrene (PS)

### Purity & certification

Satisfies the requirement IATA, ADR  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Purity standard sterile  
Sterilisation Electron irradiation  
Batched yes

### Packaging

Minimum order qty. 500  
Type of smallest subpackaging case  
Piece(s) / inner box 50  
Piece (s) / outer case 500  
Piece(s) / pallet 24000  
Depth of box 204 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® Lithium heparin LH, 9 ml, cap orange, (LxØ): 92 x 16 mm, with plastic label**

**Page 3**

Width of box	134 mm
Height of box	88 mm
Depth of case	648 mm
Width of case	303 mm
Height of case	150 mm
Case volume	0.0295 cbm
Weight of product	0.0088 kg
Weight of case	5.33 kg
EAN of inner box	4038917085236
EAN case	4038917003926

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**