



Product description

Order number	03.1631.001
Product description	S-Monovette® Lithium heparin gel LH, preparation: Lithium heparin gel, 4.7 ml, membrane screw cap, cap orange, colour code EU, (LxØ) without cap: 75 x 15 mm, with paper label, label/print: orange, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Lithium heparin gel
Preparation concentration	25 IU/ml blood
Colour code	EU
Label/ Print	with paper label
Colour of print/label	orange
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2500 x g
Centrifugation temperature	20
Centrifugation time	15 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Liquid-tight Connection	yes membrane screw cap
-------------------------	---------------------------

Size

Sample volume	4.7 ml
Outer diameter of cap	18 mm
Diameter	15 mm
Length including cap	91 mm
Length excluding cap	75 mm
Length including cap and plunger	112 mm

Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	orange
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	26000
Depth of box	189 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Lithium heparin gel LH, 4.7 ml, cap orange, (LxØ): 75 x 15 mm, with paper label

Width of box	117 mm
Height of box	82 mm
Depth of case	440 mm
Width of case	400 mm
Height of case	129 mm
Case volume	0.0227 cbm
Weight of product	0.0075 kg
Weight of case	4.54 kg
EAN of inner box	4038917132206
EAN case	4038917024310

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

