



## Product Specification

**S-Monovette® EDTA K3E, 2.6 ml, cap red, (LxØ):  
65 x 13 mm, with plastic label pre-barcoded,  
Pre-barcode with 8-digit unique number range  
and 3-digit prefix**

Page 1



### Product description

Order number	04.1901.088
Product description	S-Monovette® EDTA K3E, preparation: K3 EDTA, 2.6 ml, membrane screw cap, cap red, colour code EU, (LxØ) without cap: 65 x 13 mm, with plastic label pre-barcoded, Pre-barcode with 8-digit unique number range and 3-digit prefix 011, label/print: white/red, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	K3 EDTA
Preparation concentration	1.6 mg/ml blood
Colour code	EU
Label/ Print	with plastic label pre-barcoded
Colour of print/label	white/red
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® EDTA K3E, 2.6 ml, cap red, (LxØ):  
65 x 13 mm, with plastic label pre-barcoded,  
Pre-barcode with 8-digit unique number range  
and 3-digit prefix**

Page 2

Liquid-tight yes  
Connection membrane screw cap

### Size

Sample volume 2.6 ml  
Outer diameter of cap 16 mm  
Diameter 13 mm  
Length including cap 81 mm  
Length excluding cap 65 mm  
Length including cap and plunger 102 mm

### Material & colours

Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap red  
Membrane material rubber  
Tube material Polypropylene (PP)  
Piston material High Density Polyethylene (HD-PE)  
Piston plunger material Polystyrene (PS)

### Purity & certification

Satisfies the requirement IATA, ADR  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Purity standard sterile  
Sterilisation Electron irradiation  
Batched yes

### Packaging

Minimum order qty. 20000  
Type of smallest case  
subpackaging  
Piece(s) / inner box 50  
Piece (s) / outer case 500  
Piece(s) / pallet 42000

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® EDTA K3E, 2.6 ml, cap red, (LxØ):  
65 x 13 mm, with plastic label pre-barcoded,  
Pre-barcode with 8-digit unique number range  
and 3-digit prefix**

**Page 3**

Depth of box	171 mm
Width of box	105 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	125 mm
Case volume	0.0174 cbm
Weight of product	0.004 kg
Weight of case	2.69 kg
EAN of inner box	4038917947947
EAN case	4038917947961

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**