



Product Specification

S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 2.9 ml, cap blue, (LxØ): 65 x 13 mm, with paper label

Page 1



Product description

Order number	04.1902.100
Product description	S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, preparation: Citrate 9NC, 2.9 ml, membrane screw cap, cap blue, colour code ISO, (LxØ) without cap: 65 x 13 mm, with paper label, label/print: blue, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Citrate 9NC
Preparation concentration	0.106 mol/l
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	blue
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	1800 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-01

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 2.9 ml, cap blue, (LxØ): 65 x 13 mm, with paper label

Page 3

Width of box	105 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	125 mm
Case volume	0.0174 cbm
Weight of product	0.0043 kg
Weight of case	2.84 kg
EAN of inner box	4038917134194
EAN case	4038917025508

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-01

This document was prepared by EDP support and is valid without signature.



SARSTEDT