



Product Specification

S-Monovette® Fluoride/EDTA FE, 2.6 ml, cap yellow, (LxØ): 65 x 13 mm, with plastic label

Page 1



Product description

Order number	04.1903
Product description	S-Monovette® Fluoride/EDTA FE, preparation: Fluoride + EDTA, 2.6 ml, membrane screw cap, cap yellow, colour code EU, (LxØ) without cap: 65 x 13 mm, with plastic label, label/print: transparent/yellow, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Fluoride + EDTA
Preparation concentration	1.0 mg/ml blood + 1.2 mg/ml blood
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	transparent/yellow
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

S-Monovette® Fluoride/EDTA FE, 2.6 ml, cap yellow, (LxØ): 65 x 13 mm, with plastic label

Page 2

Liquid-tight Connection yes
membrane screw cap

Size

Sample volume 2.6 ml
Outer diameter of cap 16 mm
Diameter 13 mm
Length including cap 81 mm
Length excluding cap 65 mm
Length including cap and plunger 102 mm

Material & colours

Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap yellow
Membrane material rubber
Tube material Polypropylene (PP)
Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 500
Type of smallest subpackaging case
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 42000
Depth of box 171 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

S-Monovette® Fluoride/EDTA FE, 2.6 ml, cap yellow, (LxØ): 65 x 13 mm, with plastic label

Page 3

Width of box	105 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	125 mm
Case volume	0.0174 cbm
Weight of product	0.004 kg
Weight of case	2.69 kg
EAN of inner box	4038917087971
EAN case	4038917006132

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT