



**Product Specification**  
**S-Monovette® Sodium heparin NH, 2.6 ml, cap green, (LxØ): 65 x 13 mm, with paper label**



**Product description**

Order number	04.1913.100
Product description	S-Monovette® Sodium heparin NH, preparation: Sodium heparin, 2.6 ml, membrane screw cap, cap green, colour code ISO, (LxØ) without cap: 65 x 13 mm, with paper label, label/print: green, 50 piece(s)/case, sterile

**Product characteristics**

Type of collection	venous
Type of preparation	Sodium heparin
Preparation concentration	20 IU/ml blood
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	green
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**S-Monovette<sup>®</sup> Sodium heparin NH, 2.6 ml, cap green, (LxØ): 65 x 13 mm, with paper label**

Page 2

Liquid-tight yes  
Connection membrane screw cap

### Size

Sample volume 2.6 ml  
Outer diameter of cap 16 mm  
Diameter 13 mm  
Length including cap 81 mm  
Length excluding cap 65 mm  
Length including cap and plunger 102 mm

### Material & colours

Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap green  
Membrane material rubber  
Tube material Polypropylene (PP)  
Piston material High Density Polyethylene (HD-PE)  
Piston plunger material Polystyrene (PS)

### Purity & certification

Satisfies the requirement IATA, ADR  
Product category In vitro diagnostic device, CE  
CE certified CE  
Purity standard sterile  
Sterilisation Electron irradiation  
Batched yes

### Packaging

Minimum order qty. 500  
Type of smallest case  
subpackaging  
Piece(s) / inner box 50  
Piece (s) / outer case 500  
Piece(s) / pallet 42000  
Depth of box 171 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**S-Monovette® Sodium heparin NH, 2.6 ml, cap green, (LxØ): 65 x 13 mm, with paper label**

Width of box	105 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	125 mm
Case volume	0.0174 cbm
Weight of product	0.0043 kg
Weight of case	2.86 kg
EAN of inner box	4038917134330
EAN case	4038917072533

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

