



Product description

| | |
|---------------------|---|
| Order number | 04.1918.100 |
| Product description | S-Monovette® Fluoride/EDTA FE, preparation: Fluoride + EDTA, 2.7 ml, membrane screw cap, cap grey, colour code ISO, (LxØ) without cap: 75 x 13 mm, with paper label, label/print: grey, 50 piece(s)/case, sterile |

Product characteristics

| | |
|----------------------------|-----------------------------------|
| Type of collection | venous |
| Type of preparation | Fluoride + EDTA |
| Preparation concentration | 1.0 mg/ml blood + 1.2 mg/ml blood |
| Colour code | ISO |
| Label/ Print | with paper label |
| Colour of print/label | grey |
| Graduation | fill mark |
| Cap | membrane screw cap |
| Closure type | screw cap |
| Recommended centrifugation | 2000 x g |
| Centrifugation temperature | 20 |
| Centrifugation time | 10 Min |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





| | |
|--------------|--------------------|
| Liquid-tight | yes |
| Connection | membrane screw cap |

Size

| | |
|----------------------------------|--------|
| Sample volume | 2.7 ml |
| Outer diameter of cap | 16 mm |
| Diameter | 13 mm |
| Length including cap | 91 mm |
| Length excluding cap | 75 mm |
| Length including cap and plunger | 112 mm |

Material & colours

| | |
|-------------------------|-----------------------------------|
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | grey |
| Membrane material | rubber |
| Tube material | Polypropylene (PP) |
| Piston material | High Density Polyethylene (HD-PE) |
| Piston plunger material | Polystyrene (PS) |

Purity & certification

| | |
|---------------------------|-------------------------------|
| Satisfies the requirement | IATA, ADR |
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Purity standard | sterile |
| Sterilisation | Electron irradiation |
| Batched | yes |

Packaging

| | |
|-------------------------------|--------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 39000 |
| Depth of box | 171 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.



Product Specification

**S-Monovette® Fluoride/EDTA FE, 2.7 ml, cap grey,
(LxØ): 75 x 13 mm, with paper label**

Page 3

| | |
|-------------------|---------------|
| Width of box | 116 mm |
| Height of box | 73 mm |
| Depth of case | 388 mm |
| Width of case | 358 mm |
| Height of case | 135 mm |
| Case volume | 0.0188 cbm |
| Weight of product | 0.0043 kg |
| Weight of case | 2.88 kg |
| EAN of inner box | 4038917228039 |
| EAN case | 4038917228022 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT