



## Product Specification

**S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 3 ml, cap green, (LxØ): 75 x 13 mm, with plastic label**

**Page 1**



### Product description

Order number	04.1919
Product description	S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, preparation: Citrate 9NC, 3 ml, membrane screw cap, cap green, colour code EU, (LxØ) without cap: 75 x 13 mm, with plastic label, label/print: transparent/green, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Citrate 9NC
Preparation concentration	0.106 mol/l
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	transparent/green
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	1800 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**





## Product Specification

**S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 3 ml,  
cap green, (LxØ): 75 x 13 mm, with plastic label**

**Page 3**

Width of box	116 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.0045 kg
Weight of case	2.99 kg
EAN of inner box	4038917223706
EAN case	4038917223690

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**