



Product Specification

**S-Monovette® Serum Gel CAT, 4.9 ml, cap brown,
(LxØ): 90 x 13 mm, with plastic label
pre-barcoded, Pre-barcode with 8-digit unique
number range and 3-digit prefix**

Page 1



Product description

Order number	04.1935.088
Product description	S-Monovette® Serum Gel CAT, preparation: Clotting activator/gel, 4.9 ml, membrane screw cap, cap brown, colour code EU/ISO, (LxØ) without cap: 90 x 13 mm, with plastic label pre-barcoded, Pre-barcode with 8-digit unique number range and 3-digit prefix 009, label/print: white/brown, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Clotting activator/gel
Colour code	EU/ISO
Label/ Print	with plastic label pre-barcoded
Colour of print/label	white/brown
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**S-Monovette® Serum Gel CAT, 4.9 ml, cap brown,
(LxØ): 90 x 13 mm, with plastic label
pre-barcoded, Pre-barcode with 8-digit unique
number range and 3-digit prefix**

Page 2

Recommended centrifugation	2500 x g
Centrifugation time	10 Min
Liquid-tight	yes
Connection	membrane screw cap

Size

Sample volume	4.9 ml
Outer diameter of cap	13 mm
Diameter	13 mm
Length including cap	106 mm
Length excluding cap	90 mm
Length including cap and plunger	127 mm

Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	brown
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**S-Monovette® Serum Gel CAT, 4.9 ml, cap brown,
(LxØ): 90 x 13 mm, with plastic label
pre-barcoded, Pre-barcode with 8-digit unique
number range and 3-digit prefix**

Page 3

Packaging

Minimum order qty.	20000
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	36000
Depth of box	171 mm
Width of box	130 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	150 mm
Case volume	0.0208 cbm
Weight of product	0.0057 kg
Weight of case	3.62 kg
EAN of inner box	4038917947992
EAN case	4038917948005

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT