



Product description

Order number	04.1936.001
Product description	S-Monovette® Lithium heparin LH, preparation: Lithium heparin, 4.9 ml, membrane screw cap, cap orange, colour code EU, (LxØ) without cap: 90 x 13 mm, with paper label, label/print: orange, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Lithium heparin
Preparation concentration	16 IU/ml blood
Colour code	EU
Label/ Print	with paper label
Colour of print/label	orange
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Lithium heparin LH, 4.9 ml, cap orange, (LxØ): 90 x 13 mm, with paper label

Liquid-tight Connection yes
 membrane screw cap

Size

Sample volume 4.9 ml
 Outer diameter of cap 16 mm
 Diameter 13 mm
 Length including cap 106 mm
 Length excluding cap 90 mm
 Length including cap and plunger 127 mm

Material & colours

Colour of product transparent
 Closure material High Density Polyethylene (HD-PE)
 Colour of cap orange
 Membrane material rubber
 Tube material Polypropylene (PP)
 Piston material High Density Polyethylene (HD-PE)
 Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
 Product category In vitrodiagnostic device, CE
 CE certified CE
 Purity standard sterile
 Sterilisation Electron irradiation
 Batched yes

Packaging

Minimum order qty. 500
 Type of smallest subpackaging case
 Piece(s) / inner box 50
 Piece (s) / outer case 500
 Piece(s) / pallet 36000
 Depth of box 171 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Lithium heparin LH, 4.9 ml, cap orange, (LxØ): 90 x 13 mm, with paper label

Width of box	130 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	150 mm
Case volume	0.0208 cbm
Weight of product	0.0054 kg
Weight of case	3.45 kg
EAN of inner box	4038917134743
EAN case	4038917025911

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

