



## Product Specification

**S-Monovette® EDTA K3E, 1.8 ml, cap red, (LxØ):  
65 x 13 mm, with plastic label pre-barcoded,  
Pre-barcode with 8-digit unique number range  
and 3-digit prefix**

**Page 1**



### Product description

|                     |                                                                                                                                                                                                                                                                                               |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Order number        | 04.1951.088                                                                                                                                                                                                                                                                                   |
| Product description | S-Monovette® EDTA K3E, K3E, preparation: K3 EDTA, 1.8 ml, membrane screw cap, cap red, colour code EU, (LxØ) without cap: 65 x 13 mm, with plastic label pre-barcoded, Pre-barcode with 8-digit unique number range and 3-digit prefix 013, label/print: white/red, 50 piece(s)/case, sterile |

### Product characteristics

|                           |                                 |
|---------------------------|---------------------------------|
| Type of collection        | venous                          |
| Type of preparation       | K3 EDTA                         |
| Preparation concentration | 1.6 mg/ml blood                 |
| Colour code               | EU                              |
| Label/ Print              | with plastic label pre-barcoded |
| Colour of print/label     | white/red                       |
| Graduation                | fill mark                       |
| Cap                       | membrane screw cap              |
| Closure type              | screw cap                       |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® EDTA K3E, 1.8 ml, cap red, (LxØ):  
65 x 13 mm, with plastic label pre-barcoded,  
Pre-barcode with 8-digit unique number range  
and 3-digit prefix**

**Page 2**

Liquid-tight Connection      yes  
membrane screw cap

### Size

Sample volume      1.8 ml  
Outer diameter of cap      16 mm  
Diameter      13 mm  
Length including cap      81 mm  
Length excluding cap      65 mm  
Length including cap and plunger      102 mm

### Material & colours

Colour of product      transparent  
Closure material      High Density Polyethylene (HD-PE)  
Colour of cap      red  
Membrane material      rubber  
Tube material      Polypropylene (PP)  
Piston material      High Density Polyethylene (HD-PE)  
Piston plunger material      Polystyrene (PS)

### Purity & certification

Satisfies the requirement      IATA, ADR  
Product category      In vitrodiagnostic device, CE  
CE certified      CE  
Purity standard      sterile  
Sterilisation      Electron irradiation  
Batched      yes

### Packaging

Minimum order qty.      20000  
Type of smallest subpackaging      case  
Piece(s) / inner box      50  
Piece (s) / outer case      500  
Piece(s) / pallet      42000

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**S-Monovette® EDTA K3E, 1.8 ml, cap red, (LxØ):  
65 x 13 mm, with plastic label pre-barcoded,  
Pre-barcode with 8-digit unique number range  
and 3-digit prefix**

|                   |               |
|-------------------|---------------|
| Depth of box      | 171 mm        |
| Width of box      | 105 mm        |
| Height of box     | 73 mm         |
| Depth of case     | 388 mm        |
| Width of case     | 358 mm        |
| Height of case    | 125 mm        |
| Case volume       | 0.0174 cbm    |
| Weight of product | 0.0039 kg     |
| Weight of case    | 2.63 kg       |
| EAN of inner box  | 4038917949385 |
| EAN case          | 4038917949392 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

