



**Product Specification**  
**S-Monovette® Hirudin, 1.6 ml, cap dark green,**  
**(LxØ): 75 x 13 mm, with paper label**

**Page 1**



### Product description

Order number	04.1959.001
Product description	S-Monovette® Hirudin, For platelet function diagnostics on the ROCHE Multiplate® , preparation: Hirudin, 1.6 ml, membrane screw cap, cap dark green, (LxØ) without cap: 75 x 13 mm, with paper label, label/print: dark green, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Hirudin
Preparation concentration	525 ATU Hirudin/ ml blood
Label/ Print	with paper label
Colour of print/label	dark green
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Example of use	For platelet function diagnostics on the ROCHE Multiplate®
Liquid-tight	yes
Connection	membrane screw cap

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**S-Monovette® Hirudin, 1.6 ml, cap dark green,**  
**(LxØ): 75 x 13 mm, with paper label**

**Page 2**

**Size**

Sample volume	1.6 ml
Outer diameter of cap	16 mm
Diameter	13 mm
Length including cap	91 mm
Length excluding cap	75 mm
Length including cap and plunger	112 mm

**Material & colours**

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	dark green
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

**Purity & certification**

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Minimum order qty.	50
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	39000
Depth of box	171 mm
Width of box	116 mm
Height of box	73 mm
Depth of case	388 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**S-Monovette® Hirudin, 1.6 ml, cap dark green,**  
**(LxØ): 75 x 13 mm, with paper label**

**Page 3**

Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.0047 kg
Weight of case	3.09 kg
EAN of inner box	4038917376983
EAN case	4038917376976

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**