



### Product description

Order number	05.1076
Product description	S-Monovette® Fluoride/heparin FH, preparation: Fluoride + heparin, 2.7 ml, membrane screw cap, cap yellow, colour code EU, (LxØ) without cap: 66 x 11 mm, with plastic label, label/print: transparent/yellow, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Fluoride + heparin
Preparation concentration	1.0 mg/ml blood + 16 IU Heparin/ml blood
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	transparent/yellow
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	9 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.







## Product Specification

**S-Monovette® Fluoride/heparin FH, 2.7 ml, cap yellow, (LxØ): 66 x 11 mm, with plastic label**

**Page 3**

Width of box	108 mm
Height of box	63 mm
Depth of case	333 mm
Width of case	303 mm
Height of case	125 mm
Case volume	0.0126 cbm
Weight of product	0.0042 kg
Weight of case	2.62 kg
EAN of inner box	4038917085281
EAN case	4038917003988

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**