



Product description

Order number	05.1553
Product description	S-Monovette® Lithium heparin LH, preparation: Lithium heparin, 2.7 ml, membrane screw cap, cap orange, colour code EU, (LxØ) without cap: 66 x 11 mm, with plastic label, label/print: transparent/orange, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Lithium heparin
Preparation concentration	16 IU/ml blood
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	transparent/orange
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Lithium heparin LH, 2.7 ml, cap orange, (LxØ): 66 x 11 mm, with plastic label

Width of box	108 mm
Height of box	63 mm
Depth of case	333 mm
Width of case	303 mm
Height of case	125 mm
Case volume	0.0126 cbm
Weight of product	0.0041 kg
Weight of case	2.61 kg
EAN of inner box	4038917087018
EAN case	4038917005500

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

