



Product Specification
**S-Monovette® neutral Z, 2.7 ml, cap white, (LxØ):
66 x 11 mm, with paper label**

Page 1



Product description

Order number	05.1729.001
Product description	S-Monovette® neutral Z, without preparation, 2.7 ml, membrane screw cap, cap white, colour code ISO, (LxØ) without cap: 66 x 11 mm, with paper label, label/print: black, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	without preparation
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	black
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min
Coding colour	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**S-Monovette® neutral Z, 2.7 ml, cap white, (LxØ):
66 x 11 mm, with paper label**

Liquid-tight Connection yes
Connection membrane screw cap

Size

Sample volume 2.7 ml
Outer diameter of cap 13 mm
Diameter 11 mm
Length including cap 82 mm
Length excluding cap 66 mm
Length including cap and plunger 103 mm

Material & colours

Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap white
Membrane material rubber
Tube material Polypropylene (PP)
Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 500
Type of smallest subpackaging case
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 56000
Depth of box 145 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® neutral Z, 2.7 ml, cap white, (LxØ):
66 x 11 mm, with paper label

Page 3

Width of box	108 mm
Height of box	63 mm
Depth of case	333 mm
Width of case	303 mm
Height of case	125 mm
Case volume	0.0126 cbm
Weight of product	0.0038 kg
Weight of case	2.45 kg
EAN of inner box	4038917133371
EAN case	4038917024853

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT