



Product Specification
Monovette® Luer Serum CAT, 4.5 ml, cap red,
(LxØ): 92 x 11 mm, with paper label

Page 1



Product description

Order number	05.274.160
Product description	Monovette® Luer Serum CAT, preparation: Clotting activator, nominal volume: 4.5 ml, (LxØ): 92 x 11 mm, Luer screw cap, cap red, with paper label, label/print: red, animal venous blood collection, 100 piece(s)/bag

Product characteristics

Type of collection	venous
Type of preparation	Clotting activator
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	red
Cap	Luer screw cap
Closure type	screw cap
Application	animal venous blood collection
Recommended centrifugation	2000 x g
Centrifugation time	10 Min
Connection	Luer

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Monovette® Luer Serum CAT, 4.5 ml, cap red,
(LxØ): 92 x 11 mm, with paper label

Page 2

Size

Sample volume	4.5 ml
Diameter	11 mm
Height of product	92 mm
Length of product	92 mm
Length including cap	117 mm
Length excluding cap	92 mm
Length including cap and plunger	127 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	red
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	100
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	500
Piece(s) / pallet	36000
Case volume	0.0182 cbm
Weight of product	0.0067 kg
Weight of case	3.6 kg
EAN of inner box	4038917104753
EAN case	4038917032322

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

