



Product Specification
S-Monovette® EDTA K3E, 1.2 ml, cap red, (LxØ):
66 x 8 mm, with plastic label

Page 1



Product description

Order number	06.1662.001
Product description	S-Monovette® EDTA K3E, preparation: K3 EDTA, 1.2 ml, without mixing ball, membrane screw cap, cap red, colour code EU, (LxØ) without cap: 66 x 8 mm, with plastic label, label/print: white/red, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	K3 EDTA
Preparation concentration	1.6 mg/ml blood
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	white/red
Graduation	yes
Cap	membrane screw cap
Closure type	screw cap
Liquid-tight	yes
Connection	membrane screw cap

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
S-Monovette® EDTA K3E, 1.2 ml, cap red, (LxØ):
66 x 8 mm, with plastic label

Page 2

Size

Sample volume	1.2 ml
Outer diameter of cap	13 mm
Diameter	8 mm
Length including cap	82 mm
Length excluding cap	66 mm
Length including cap and plunger	103 mm

Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	red
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	rubber
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	56000
Depth of box	145 mm
Width of box	108 mm
Height of box	63 mm
Depth of case	333 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
S-Monovette® EDTA K3E, 1.2 ml, cap red, (LxØ):
66 x 8 mm, with plastic label

Width of case	303 mm
Height of case	125 mm
Case volume	0.0126 cbm
Weight of product	0.0031 kg
Weight of case	2.12 kg
EAN of inner box	4038917259828
EAN case	4038917259811

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

