



Product description

Order number	06.1669.001
Product description	S-Monovette® Lithium heparin gel LH, preparation: Lithium heparin gel, 1.1 ml, membrane screw cap, cap orange, colour code EU, (LxØ) without cap: 66 x 8 mm, with plastic label, label/print: white/orange, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Lithium heparin gel
Preparation concentration	35 IU/ml blood
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	white/orange
Graduation	yes
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2500 x g
Centrifugation temperature	20
Centrifugation time	15 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Lithium heparin gel LH, 1.1 ml, cap orange, (LxØ): 66 x 8 mm, with plastic label

Liquid-tight yes
Connection membrane screw cap

Size

Sample volume 1.1 ml
Outer diameter of cap 13 mm
Diameter 8 mm
Length including cap 82 mm
Length excluding cap 66 mm
Length including cap and plunger 103 mm

Material & colours

Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap orange
Membrane material rubber
Tube material Polypropylene (PP)
Piston material rubber
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 500
Type of smallest case
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 56000
Depth of box 145 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Lithium heparin gel LH, 1.1 ml, cap orange, (LxØ): 66 x 8 mm, with plastic label

Page 3

Width of box	108 mm
Height of box	63 mm
Depth of case	333 mm
Width of case	303 mm
Height of case	125 mm
Case volume	0.0126 cbm
Weight of product	0.0033 kg
Weight of case	2.18 kg
EAN of inner box	4038917132930
EAN case	4038917079396

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT