



### **Product description**

Order number	74.5330.500
Product description	UD 3plus, Urine drainage system, 2 l, with holder, sterile, tube length: 100 cm, 1 piece(s)/blister
PZN number	16086848

### **Product characteristics**

Colour of print/label	blue
Design	with holder
Application	Urine drainage
Volume of work	2 l

### **Size**

Sample volume	2 l
Max. volume	2 l
Diameter of tube	10.1 mm
Length of tube	100 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material Polyvinylchloride (PVC)  
Colour of product transparent

**Purity & certification**

Product category Medical device  
CE certified CE supplier  
Purity standard sterile  
Sterilisation Gas sterilisation  
Batched yes

**Packaging**

Minimum order qty. 20  
Type of smallest blister  
subpackaging  
Piece(s) / inner box 1  
Piece (s) / outer case 20  
Piece(s) / pallet 1120  
Depth of case 395 mm  
Width of case 295 mm  
Height of case 220 mm  
Case volume 0.0256 cbm  
Weight of product 0.191 kg  
Weight of case 3.82 kg  
EAN of inner box 6955824024333  
EAN case 6955824024180

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

