



Product description

Order number	74.5684.002
Product description	Urine measuring system exchange bag, 2 l, with drainage valve, sterile, 1 piece(s)/blister

Product characteristics

Colour of print/label	blue
Design	with drainage valve
Application	Urine drainage
Volume of work	2 l

Size

Sample volume	2 l
Max. volume	2 l

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polyvinylchloride (PVC)
Colour of product transparent

Purity & certification

Product category Medical device
CE certified CE supplier
Purity standard sterile
Sterilisation Gas sterilisation
Batched yes

Packaging

Minimum order qty. 50
Type of smallest blister
subpackaging
Piece(s) / inner box 1
Piece (s) / outer case 50
Piece(s) / pallet 4050
Depth of case 395 mm
Width of case 245 mm
Height of case 200 mm
Case volume 0.0194 cbm
Weight of product 0.0401 kg
Weight of case 2.77 kg
EAN of inner box 6955824024449
EAN case 6955824024067

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

