




# Product Specification

## Urine bag, 1.5 l, with anti-reflux valve, with drainage valve

Page 1



### Product description

Order number	74.5218.101
Product description	Urine bag, 1.5 l, with anti-reflux valve, with drainage valve, tube length: 90 cm, 150 piece(s)/case
PZN number	16333689

### Product characteristics

Colour of print/label	blue
Design	with drainage valve
Application	Urine drainage
Volume of work	1,5 l
Anti-reflux valve	yes

### Size

Sample volume	1.5 l
Max. volume	1.5 l
Diameter of tube	6.6 mm
Length of tube	90 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Urine bag, 1.5 l, with anti-reflux valve, with drainage valve**

**Page 2**

**Material & colours**

Product material Polyvinylchloride (PVC)  
Colour of product transparent

**Purity & certification**

Product category Medical device  
CE certified CE supplier  
Batched yes

**Packaging**

Minimum order qty. 150  
Type of smallest case  
subpackaging  
Piece(s) / inner box 150  
Piece (s) / outer case 150  
Piece(s) / pallet 9600  
Depth of case 580 mm  
Width of case 390 mm  
Height of case 150 mm  
Case volume 0.0339 cbm  
Weight of product 0.0352 kg  
Weight of case 5.8 kg  
EAN of inner box 6971243660969  
EAN case 6971243660693

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**