



### **Product description**

Order number	79.8710.577
Product description	Ink ribbon, Consumables for protocol printer SP742MD, black

### **Product characteristics**

Application	Warming
Suitable for	Protocol printer

### **Size**

Width of product	7 cm
Height of product	3 cm
Length of product	12 cm

### **Material & colours**

Colour of product	black
-------------------	-------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

### **Packaging**

Minimum order qty.                      1  
Type of smallest                              case  
subpackaging  
Piece(s) / inner box                      1  
Piece (s) / outer case                      1  
Piece(s) / pallet                              32  
Depth of case                              440 mm  
Width of case                                390 mm  
Height of case                               220 mm  
Case volume                                0.0378 cbm  
Weight of product                        3.7 kg  
Weight of case                              4.5 kg  
EAN of inner box                        4038917975940  
EAN case                                    4038917975957

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

