



Product description

Order number 97.8710.541
Product description Service module SAHARA, for adjusting and calibrating the SAHARA-III, TRANSTHERM-III and SAHARA-TSC

Product characteristics

Application Warming

Size

Width of product 28 cm
Height of product 7.5 cm
Length of product 33 cm

Purity & certification

Product category no medical device | IVD
Batched no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1
Case volume	0.0314 cbm
Weight of product	1.85 kg
Weight of case	2.18 kg
EAN of inner box	4028339753226
EAN case	4028339753219

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

