



Product description

| | |
|---------------------|--|
| Order number | 94.6001.153 |
| Product description | miniPERM® , Nutrient module, for miniPERM® production module, autoclavable up to 121 °C, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag |

Product characteristics

| | |
|------------------------|-----------------------------|
| Closure type | screw cap |
| Suitable for | miniPERM® production module |
| Volume of work | 350 ml - 450 ml |
| Autoclavable (max. °C) | 121 °C |

Size

| | |
|-------------------|--------|
| Sample volume | 400 ml |
| Diameter | 102 mm |
| Height of product | 129 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|------------------|-----------------------------------|
| Product material | Polycarbonate (PC) |
| Closure material | High Density Polyethylene (HD-PE) |

Purity & certification

| | |
|------------------|---|
| Product category | no medical device IVD |
| Purity standard | non-pyrogenic/endotoxin-free, non-cytotoxic |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 4 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 4 |
| Piece (s) / outer case | 24 |
| Piece(s) / pallet | 480 |
| Depth of case | 498 mm |
| Width of case | 364 mm |
| Height of case | 310 mm |
| Case volume | 0.0562 cbm |
| Weight of product | 0.1208 kg |
| Weight of case | 4.74 kg |
| EAN of inner box | 4038917254991 |
| EAN case | 4038917254984 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

