



### Product description

Order number	83.3944.001
Product description	Pressure filtration unit, Filtropur L 0.2 LS, pre-filter: GF, membrane: PES, filtration surface: 20 cm <sup>2</sup> , pore size: 0.2 µm, for sterile filtration, connection: Luer lock, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag

### Product characteristics

Example of use	for sterile filtration
Filter	yes
Surface area of filter	20 cm <sup>2</sup>
Pore size	0,2
Connection	Luer lock

### Size

Filtration volume	10000 ml
-------------------	----------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polycarbonate (PC)
Membrane material	Polyethersulfone (PES)

### **Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gas sterilisation
Non-cytotoxic	yes
Batched	yes

### **Packaging**

Minimum order qty.	10
Type of smallest subpackaging	bag
Piece(s) / inner box	10
Piece (s) / outer case	50
Piece(s) / pallet	400
Depth of box	352 mm
Width of box	155 mm
Height of box	123 mm
Depth of case	530 mm
Width of case	370 mm
Height of case	280 mm
Case volume	0.0549 cbm
Weight of product	0.025 kg
Weight of case	4.24 kg
EAN of inner box	4038917926188
EAN case	4038917926195

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

