



Product Specification
**CryoPure tubes, 2 ml, QuickSeal screw cap,
yellow**

Page 1



Product description

| | |
|---------------------|--|
| Order number | 72.379.004 |
| Product description | CryoPure tubes, 2 ml, tube: PP, QuickSeal screw cap, cap assembled, HD-PE, yellow, external thread, nominal volume: 1.8 ml, Cryo Performance Tested, 50 piece(s)/bag |

Product characteristics

| | |
|---------------------------------------|---------------------|
| Label/ Print | with print |
| Colour of print/label | white |
| Graduation | yes |
| Cap | QuickSeal screw cap |
| Closure type | screw cap |
| Internally or externally threaded cap | external thread |
| Volume of work | 1,8 ml |
| Vessel type | Tube |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**CryoPure tubes, 2 ml, QuickSeal screw cap,
yellow**

Page 2

Size

| | |
|----------------------|----------|
| Sample volume | 1.8 ml |
| Max. volume | 2 ml |
| Diameter | 12.65 mm |
| Length of product | 46.7 mm |
| Length excluding cap | 39.8 mm |

Material & colours

| | |
|------------------|-----------------------------------|
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | yellow |
| Tube material | Polypropylene (PP) |

Purity & certification

| | |
|------------------|--------------------------------|
| Product category | In vitro diagnostic device, CE |
| CE certified | CE |
| Purity standard | Cryo Performance Tested |
| Sterilisation | Electron irradiation |
| Pyrogen-free | yes |
| Non-cytotoxic | yes |
| Batched | yes |

Packaging

| | |
|-------------------------------|------------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 2000 |
| Piece(s) / pallet | 104000 |
| Depth of box | 293 mm |
| Width of box | 183 mm |
| Height of box | 118 mm |
| Depth of case | 598 mm |
| Width of case | 378 mm |
| Height of case | 135 mm |
| Case volume | 0.0305 cbm |
| Weight of product | 0.002 kg |
| Weight of case | 5.17 kg |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**CryoPure tubes, 2 ml, QuickSeal screw cap,
yellow**

Page 3

| | |
|------------------|---------------|
| EAN of inner box | 4038917106047 |
| EAN case | 4038917080576 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-06-02

This document was prepared by EDP support and is valid without signature.



SARSTEDT