



Product description

Order number 93.872.425
Product description Cryobox, with numerical coding at each aperture, for low-temperature storage, material: PC, yellow, slip-on lid with ventilation function, cap: transparent, (LxWxH): 75 x 75 x 52 mm, format: 5 x 5, for 25 collection tubes, for CryoPure tubes 1.2?2.0 ml internal and external thread, 5 piece(s)/bag

Product characteristics

Graduation no
Cap slip-on lid with ventilation function
Closure type slip-on lid
Design with numerical coding at each aperture
Suitable for CryoPure tubes 1.2?2.0 ml internal and external thread
Storage capacity 25
Format 5 x 5

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

Width of product	75 mm
Height of product	52 mm
Length of product	75 mm

Material & colours

Product material	Polycarbonate (PC)
Colour of product	yellow
Colour of cap	transparent

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	5
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	20
Piece(s) / pallet	2640
Depth of case	156 mm
Width of case	143 mm
Height of case	502 mm
Case volume	0.0112 cbm
Weight of product	0.058 kg
Weight of case	1.41 kg
EAN of inner box	4038917122818
EAN case	4038917420174

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

