



Product description

Order number	83.3903
Product description	Tissue culture dish, (ØxH): 150 x 20 mm, material: PS, surface: Standard, for Adherent cells, coding colour: red, TC Tested, 5 piece(s)/bag

Product characteristics

Closure type	slip-on lid
Cell type	Adherent cells
Base shape	SUREGrip Design
Growth area	151.4 cm ²
Volume of work	36 ml
Vessel type	Dish
Coding colour	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tissue culture dish, (ØxH): 150 x 20 mm, surface:
Standard

Page 2

Size

Diameter	150 mm
Height of product	20 mm

Material & colours

Product material	Polystyrene (PS)
Closure material	Polystyrene (PS)

Purity & certification

Product category	no medical device IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	100
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	100
Piece(s) / pallet	1800
Depth of case	438 mm
Width of case	324 mm
Height of case	540 mm
Case volume	0.0766 cbm
Weight of product	0.0541 kg
Weight of case	7.19 kg
EAN of inner box	4038917338288
EAN case	4038917338271

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT