



### Product description

|                     |  |
|---------------------|--|
| Order number        | 83.3921  |
| Product description | Cell culture plate, 12 well, material: PS, surface: Standard, for Adherent cells, coding colour: red, flat base, TC Tested, 1 piece(s)/blister |

### Product characteristics

|                 |                      |
|-----------------|----------------------|
| Label/ Print    | with plastic label   |
| Closure type    | slip-on lid          |
| Cell type       | Adherent cells       |
| Base shape      | flat base            |
| Number of wells | 12                   |
| Growth area     | 3.65 cm <sup>2</sup> |
| Volume of work  | 2 ml                 |
| Vessel type     | Plate                |
| Coding colour   | red                  |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### Size

|                   |          |
|-------------------|----------|
| Max. volume       | 6.64 ml  |
| Width of product  | 85.2 mm  |
| Height of product | 22.5 mm  |
| Length of product | 127.8 mm |

### Material & colours

|                   |                  |
|-------------------|------------------|
| Product material  | Polystyrene (PS) |
| Colour of product | transparent      |
| Closure material  | Polystyrene (PS) |

### Purity & certification

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Purity standard  | TC Tested               |
| Sterilisation    | Gamma irradiation       |
| Pyrogen-free     | yes                     |
| Non-cytotoxic    | yes                     |
| Batched          | yes                     |

### Packaging

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 50            |
| Type of smallest subpackaging | blister       |
| Piece(s) / inner box          | 1             |
| Piece (s) / outer case        | 50            |
| Piece(s) / pallet             | 2400          |
| Depth of case                 | 332 mm        |
| Width of case                 | 297 mm        |
| Height of case                | 275 mm        |
| Case volume                   | 0.0271 cbm    |
| Weight of product             | 0.0501 kg     |
| Weight of case                | 3.25 kg       |
| EAN of inner box              | 4038917338431 |
| EAN case                      | 4038917338424 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

