



Product description

| | |
|---------------------|--|
| Order number | 83.3932.041 |
| Product description | TC insert, for 24-well plates, membrane: PET, transparent, pore size: 0.4 μm , sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister |

Product characteristics

| | |
|----------------|-------------------|
| Suitable for | 24-well plates |
| Base shape | flat base |
| Growth area | 0.3 cm^2 |
| Volume of work | 0,1 ml - 0,4 ml |
| Vessel type | TC insert |
| Pore size | 0,4 |

Size

| | |
|-------------------|----------|
| Height of product | 16.32 mm |
|-------------------|----------|

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|-----------------------------------|
| Membrane material | Polyethylene therephthalate (PET) |
| Frame material | Polystyrene (PS) |

Purity & certification

| | |
|------------------|--|
| Product category | no medical device IVD |
| Purity standard | sterile, non-pyrogenic/endotoxin-free, non-cytotoxic |
| Sterilisation | X-Ray |
| Pyrogen-free | yes |
| Non-cytotoxic | yes |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 48 |
| Type of smallest subpackaging | blister |
| Piece(s) / inner box | 48 |
| Piece (s) / outer case | 1440 |
| Piece(s) / pallet | 23040 |
| Depth of box | 196 mm |
| Width of box | 104 mm |
| Height of box | 131 mm |
| Depth of case | 410 mm |
| Width of case | 435 mm |
| Height of case | 550 mm |
| Case volume | 0.0981 cbm |
| Weight of product | 0.0032 kg |
| Weight of case | 7.87 kg |
| EAN of inner box | 4038917371346 |
| EAN case | 4038917371339 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

