



### **Product description**

Order number	83.3950
Product description	Cell scraper, 2-position blade, size: S, handle length: 24 cm, blade length: 1.35 cm, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

### **Size**

Length of product	24.7 cm
Length of blade	1.35 cm
Length of stem	24 cm
Width of blade	1.35 cm

### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	red/opaque
Stem material	Polystyrene (PS)
Blade material	Thermoplastic elastomers (TPE)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

### **Packaging**

Minimum order qty.	100
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	100
Piece(s) / pallet	16500
Depth of case	276 mm
Width of case	110 mm
Height of case	333 mm
Case volume	0.0101 cbm
Weight of product	0.0081 kg
Weight of case	1.3 kg
EAN of inner box	4038917426633
EAN case	4038917426626

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

