



### **Product description**

Order number	82.1194.500
Product description	Petri dish, (ØxH): 60 x 15 mm, material: PS, transparent, with ventilation cams, sterile, 20 piece(s)/tube bag

### **Product characteristics**

Minimum temperature	0 °C
Maximum temperature	80 °C
Ventilation cams	yes

### **Size**

Max. volume	30 ml
Diameter	60 mm
Height of product	15 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Petri dish, 60 x 15 mm, transparent, with ventilation cams**

**Page 2**

**Material & colours**

Product material Polystyrene (PS)  
Colour of product transparent

**Purity & certification**

Product category no medical device | IVD  
Latex-free yes  
Purity standard sterile  
Sterilisation Gamma irradiation  
Safety data sheet no  
Batched yes

**Packaging**

Minimum order qty. 500  
Type of smallest subpackaging tube bag  
Piece(s) / inner box 20  
Piece (s) / outer case 500  
Piece(s) / pallet 12000  
Depth of case 317 mm  
Width of case 317 mm  
Height of case 365 mm  
Case volume 0.0367 cbm  
Weight of product 0.0074 kg  
Weight of case 4.49 kg  
EAN of inner box 4038917347891  
EAN case 4038917347884

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**