




Product Specification
Petri dish, 100 x 100 x 20 mm, transparent,
without ventilation cams

Page 1



Product description

| | |
|---------------------|--|
| Order number | 82.9923.422 |
| Product description | Petri dish, (LxWxH): 100 x 100 x 20 mm, material: PS, transparent, without ventilation cams, sterile, 4 piece(s)/bag |

Product characteristics

| | |
|---------------------|-------|
| Minimum temperature | 0 °C |
| Maximum temperature | 80 °C |
| Ventilation cams | no |

Size

| | |
|-------------------|--------|
| Max. volume | 175 ml |
| Width of product | 100 mm |
| Height of product | 20 mm |
| Length of product | 100 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Petri dish, 100 x 100 x 20 mm, transparent,
without ventilation cams

Page 2

Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Latex-free yes
Purity standard sterile
Sterilisation Gamma irradiation
Safety data sheet no
Batched yes

Packaging

Minimum order qty. 160
Type of smallest
subpackaging bag
Piece(s) / inner box 4
Piece (s) / outer case 160
Piece(s) / pallet 4800
Depth of case 608 mm
Width of case 324 mm
Height of case 264 mm
Case volume 0.052 cbm
Weight of product 0.0343 kg
Weight of case 6.35 kg
EAN of inner box 4038917154178
EAN case 4038917027854

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-05

This document was prepared by EDP support and is valid without signature.



SARSTEDT