



### **Product description**

Order number	67.746
Product description	Semi-micro cuvette, max. working volume: 3 ml, material: PS, transparent, (HxW): 45 x 12 mm, optical sides: 2, path length: 10 mm, 500 piece(s)/bag

### **Product characteristics**

Volume of work	1 ml - 3 ml
Number of optical sides	2
Path length	10 mm
Wavelength in the cuvette area	330 nm

### **Size**

Max. volume	3 ml
Width of product	12 mm
Height of product	45 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Semi-micro cuvette, 3 ml, (HxW): 45 x 12 mm, PS,**  
**transparent, optical sides: 2**

**Page 2**

**Material & colours**

Product material                      Polystyrene (PS)  
Colour of product                      transparent

**Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

**Packaging**

Minimum order qty.                      2000  
Type of smallest  
subpackaging                              bag  
Piece(s) / inner box                      500  
Piece (s) / outer case                      2000  
Piece(s) / pallet                              96000  
Depth of case                              388 mm  
Width of case                                378 mm  
Height of case                               215 mm  
Case volume                                0.0315 cbm  
Weight of product                        0.0023 kg  
Weight of case                              5.02 kg  
EAN of inner box                        4038917083928  
EAN case                                    4038917002806

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**