



### **Product description**

Order number	86.1567.051
Product description	Inoculation loop, volume: 1 µl, material: PS, white, sterile, 48 piece(s)/blister

### **Size**

Width of product	4 mm
Length of product	198 mm

### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	white

### **Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Packaging**

Minimum order qty.	192000
Type of smallest subpackaging	blister
Piece(s) / inner box	960
Piece (s) / outer case	1920
Piece(s) / pallet	253440
Depth of box	216 mm
Width of box	161 mm
Height of box	101 mm
Depth of case	338 mm
Width of case	228 mm
Height of case	120 mm
Case volume	0.0092 cbm
Weight of product	0.0011 kg
Weight of case	2.53 kg
EAN of inner box	4038917258166
EAN case	4038917258159

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

