



Product description

Order number 86.1568.050
Product description Inoculation needle, material: PS, orange, sterile, 50 piece(s)/blister

Size

Width of product 4 mm
Length of product 197 mm

Material & colours

Product material Polystyrene (PS)
Colour of product orange

Purity & certification

Product category no medical device | IVD
Purity standard sterile
Sterilisation Gamma irradiation
Batched yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 2000 |
| Type of smallest subpackaging | blister |
| Piece(s) / inner box | 1000 |
| Piece (s) / outer case | 2000 |
| Piece(s) / pallet | 264000 |
| Depth of box | 216 mm |
| Width of box | 161 mm |
| Height of box | 101 mm |
| Depth of case | 338 mm |
| Width of case | 228 mm |
| Height of case | 120 mm |
| Case volume | 0.0092 cbm |
| Weight of product | 0.0012 kg |
| Weight of case | 2.79 kg |
| EAN of inner box | 4038917130974 |
| EAN case | 4038917023542 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

