



Product description

| | |
|---------------------|--|
| Order number | 93.1646.080 |
| Product description | Silicone strip, for DishRack 80, length: 370 mm, material: silicone, translucent, 1 piece(s)/bag |

Product characteristics

| | |
|--------------|-------------|
| Suitable for | DishRack 80 |
|--------------|-------------|

Size

| | |
|-------------------|--------|
| Width of product | 65 mm |
| Length of product | 370 mm |

Material & colours

| | |
|-------------------|-------------|
| Product material | silicone |
| Colour of product | translucent |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of smallest subpackaging bag
Piece(s) / inner box 1
Piece (s) / outer case 30
Piece(s) / pallet 3240
Depth of case 383 mm
Width of case 193 mm
Height of case 190 mm
Case volume 0.014 cbm
Weight of product 0.0579 kg
Weight of case 2 kg
EAN of inner box 4038917132596
EAN case 4038917040655

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

