



Product description

Order number	93.1647.050
Product description	Silicone strip, for DishRack 50, length: 240 mm, material: silicone, translucent, 1 piece(s)/bag

Product characteristics

Suitable for	DishRack 50
--------------	-------------

Size

Width of product	65 mm
Length of product	240 mm

Material & colours

Product material	silicone
Colour of product	translucent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of smallest bag
subpackaging
Piece(s) / inner box 1
Piece (s) / outer case 30
Piece(s) / pallet 4860
Depth of case 253 mm
Width of case 193 mm
Height of case 190 mm
Case volume 0.0093 cbm
Weight of product 0.0606 kg
Weight of case 2 kg
EAN of inner box 4038917132640
EAN case 4038917040648

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

