



### **Product description**

|                     |  |
|---------------------|--|
| Order number        | 80.1366  |
| Product description | Swabs, with tube included, Amies gel with charcoal, (LxØ): 175 x 12 mm, stem material: aluminium, stem length: 133 mm, swab material: viscose, sterile, 1 piece(s)/blister |

### **Product characteristics**

|                |                         |
|----------------|-------------------------|
| Closure type   | push cap                |
| Design         | with tube included      |
| Type of medium | Amies gel with charcoal |

### **Size**

|                                  |        |
|----------------------------------|--------|
| Diameter                         | 12 mm  |
| Length of product                | 175 mm |
| Length including cap             | 175 mm |
| Length excluding cap             | 155 mm |
| Length including cap and plunger | 164 mm |
| Length of stem                   | 133 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

|                  |                                   |
|------------------|-----------------------------------|
| Product material | Polypropylene (PP)                |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap    | orange                            |
| Swab material    | viscose                           |
| Stem material    | aluminium                         |

**Purity & certification**

|                   |                   |
|-------------------|-------------------|
| Product category  | Medical device    |
| CE certified      | CE supplier       |
| Purity standard   | sterile           |
| Sterilisation     | Gamma irradiation |
| Safety data sheet | yes               |
| Batched           | yes               |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 500           |
| Type of smallest subpackaging | blister       |
| Piece(s) / inner box          | 50            |
| Piece (s) / outer case        | 500           |
| Piece(s) / pallet             | 11500         |
| Depth of case                 | 520 mm        |
| Width of case                 | 450 mm        |
| Height of case                | 220 mm        |
| Case volume                   | 0.0515 cbm    |
| Weight of product             | 0.013 kg      |
| Weight of case                | 7.28 kg       |
| EAN of inner box              | 4038917086356 |
| EAN case                      | 4038917004916 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

