



Product description

| | |
|---------------------|--|
| Order number | 92.180.614 |
| Product description | Block rotor, for 7 tubes up to 28 mm diameter, for SARMIX® M 2000, 1 piece(s)/case |

Product characteristics

| | |
|------------------|----------------|
| Suitable for | SARMIX® M 2000 |
| Storage capacity | 7 |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|-------------------------------|------|
| Minimum order qty. | 1 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 10 |
| Piece(s) / pallet | 10 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Block rotor, for 7 tubes up to 28 mm diameter, for
SARMIX® M 2000

| | |
|-------------------|---------------|
| Case volume | 0.0122 cbm |
| Weight of product | 0.2009 kg |
| Weight of case | 2.21 kg |
| EAN of inner box | 4038917102025 |
| EAN case | 4038917075787 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

