



### **Product description**

Order number 92.181.611  
Product description Rotor, for 22 collection tubes

### **Product characteristics**

Storage capacity 22

### **Purity & certification**

Product category no medical device | IVD  
Batched yes

### **Packaging**

Minimum order qty. 1  
Piece(s) / inner box 1  
Piece (s) / outer case 1  
Piece(s) / pallet 40  
Case volume 0.0052 cbm  
Weight of product 1.3595 kg  
Weight of case 1.36 kg  
EAN of inner box 4038917156257  
EAN case 4038917156240

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

