




Product Specification

Closure mat, for Deep Well MegaBlock® 0.5 ml (82.1969.002) & 1.2 ml (82.1970.002 & 82.1971.002), (LxW): 79 x 121 mm, TPE, transparent

Page 1



Product description

Order number	95.1990.002
Product description	Closure mat, for Deep Well MegaBlock® 0.5 ml (82.1969.002) & 1.2 ml (82.1970.002 & 82.1971.002), (LxW): 79 x 121 mm, material: TPE, transparent, pierceable, Free from DNA/DNase/RNase, free from pyrogens/endotoxins, 10 piece(s)/bag

Product characteristics

Suitable for	Deep Well MegaBlock® 0.5 ml (82.1969.002) & 1.2 ml (82.1970.002 & 82.1971.002)
--------------	--

Size

Width of product	121 mm
Length of product	79 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

Closure mat, for Deep Well MegaBlock® 0.5 ml
(82.1969.002) & 1.2 ml (82.1970.002 &
82.1971.002), (LxW): 79 x 121 mm, TPE,
transparent

Page 2

Material & colours

Product material Thermoplastic elastomers (TPE)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Purity standard Free from DNA/DNase/RNase, free from
pyrogens/endotoxins
Batched yes

Packaging

Minimum order qty. 50
Type of smallest bag
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 250
Piece(s) / pallet 5500
Depth of case 380 mm
Width of case 280 mm
Height of case 200 mm
Case volume 0.0213 cbm
Weight of product 0.0144 kg
Weight of case 4.1 kg
EAN of inner box 4038917135023
EAN case 4038917042048

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT