



### **Product description**

Order number	82.1581.220
Product description	ELISA plate, 96 well, flat base, without lid, material: PS, black, High Binding, free from DNA/DNase/RNase, free from pyrogens and endotoxins, non-cytotoxic, 25 piece(s)/bag

### **Product characteristics**

Label/ Print	with print
Base shape	flat base
Number of wells	96
Vessel type	Plate

### **Size**

Max. volume	390 $\mu$ l
Width of product	85.2 mm
Height of product	14.35 mm
Length of product	127.8 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material Polystyrene (PS)  
Colour of product black

### **Purity & certification**

Product category no medical device | IVD  
Purity standard free from DNA/DNase/RNase, free from pyrogens and endotoxins, non-cytotoxic  
Pyrogen-free yes  
Batched yes

### **Packaging**

Minimum order qty. 50  
Type of smallest subpackaging bag  
Piece(s) / inner box 25  
Piece (s) / outer case 50  
Piece(s) / pallet 2000  
Depth of case 279 mm  
Width of case 232 mm  
Height of case 365 mm  
Case volume 0.0236 cbm  
Weight of product 0.0408 kg  
Weight of case 2.86 kg  
EAN of inner box 4038917365390  
EAN case 4038917365383

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

