



**Product Specification**  
**Micro test plate, 96 well, round base, PS,  
transparent**

**Page 1**



**Product description**

Order number	82.1582
Product description	Micro test plate, 96 well, round base, without lid, material: PS, transparent, free from DNA/DNase/RNase, free from pyrogens and endotoxins, non-cytotoxic, 25 piece(s)/bag

**Product characteristics**

Base shape	round base
Number of wells	96
Vessel type	Plate

**Size**

Max. volume	310 $\mu$ l
Width of product	85.2 mm
Height of product	14.35 mm
Length of product	127.8 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Micro test plate, 96 well, round base, PS,  
transparent**

**Page 2**

**Material & colours**

Product material Polystyrene (PS)  
Colour of product transparent

**Purity & certification**

Product category no medical device | IVD  
Purity standard free from DNA/DNase/RNase, free from pyrogens and  
endotoxins, non-cytotoxic  
Batched yes

**Packaging**

Minimum order qty. 100  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 25  
Piece (s) / outer case 100  
Piece(s) / pallet 3200  
Depth of case 359 mm  
Width of case 275 mm  
Height of case 365 mm  
Case volume 0.036 cbm  
Weight of product 0.0343 kg  
Weight of case 4.33 kg  
EAN of inner box 4038917087100  
EAN case 4038917355469

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**