



### **Product description**

Order number	83.9923.972
Product description	Micro test plate, Terasaki, 60 well, slip-on lid, base shape: conical, lid assembled, material: PS, transparent, sterile, 10 piece(s)/bag

### **Product characteristics**

Closure type	slip-on lid
Base shape	conical
Number of wells	60
Vessel type	Plate
Type of surface treatment	TC-treated

### **Size**

Max. volume	10 $\mu$ l
Width of product	57.8 mm
Height of product	9.7 mm
Length of product	85.1 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Micro test plate, Terasaki, 60 well, slip-on lid,**  
**base shape: conical, PS, transparent**

**Material & colours**

Product material Polystyrene (PS)  
Colour of product transparent

**Purity & certification**

Product category no medical device | IVD  
Purity standard sterile  
Sterilisation Gamma irradiation  
Batched yes

**Packaging**

Minimum order qty. 650  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 10  
Piece (s) / outer case 650  
Piece(s) / pallet 16250  
Depth of case 608 mm  
Width of case 324 mm  
Height of case 264 mm  
Case volume 0.052 cbm  
Weight of product 0.0173 kg  
Weight of case 12.2 kg  
EAN of inner box 4038917154710  
EAN case 4038917354363

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

