



**Product Specification**  
**Weigh tray, 250 ml, (LxW): 128 x 128 mm, PVC,  
white**

**Page 1**



**Product description**

Order number 71.9923.210  
Product description Weigh tray, 250 ml, (LxW): 128 x 128 mm, material: PVC,  
white, 100 piece(s)/bag

**Size**

Sample volume 250 ml  
Max. volume 300 ml  
Width of product 128 mm  
Length of product 128 mm

**Material & colours**

Product material Polyvinylchloride (PVC)  
Colour of product white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category no medical device | IVD  
Batched yes

**Packaging**

Minimum order qty. 1000  
Type of smallest subpackaging bag  
Piece(s) / inner box 100  
Piece (s) / outer case 1000  
Piece(s) / pallet 28000  
Depth of case 458 mm  
Width of case 373 mm  
Height of case 238 mm  
Case volume 0.0407 cbm  
Weight of product 0.005 kg  
Weight of case 5.6 kg  
EAN of inner box 4038917153645  
EAN case 4038917038225

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

