



### Product description

Order number	84.2620.002
Product description	Adhesive label, (LxW): 42 x 12 mm, material: paper, white, 4,000 piece(s)/roll

### Size

Width of product	12 mm
Length of product	42 mm

### Material & colours

Product material	paper
Colour of product	white

### Purity & certification

Product category	no medical device   IVD
Batched	yes

### Packaging

Minimum order qty.	4000
Type of smallest subpackaging	roll
Piece(s) / inner box	4000
Piece (s) / outer case	8000
Piece(s) / pallet	2592000
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.0001 kg
Weight of case	1.2 kg
EAN of inner box	4038917416917
EAN case	4038917416900

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

