



### **Product description**

|                     |   |
|---------------------|---|
| Order number        | 65.1998.400   |
| Product description | PCR strip of lids, for PCR plates (REF 72.1979.102/.132, 72.1982.202), PCR Performance Tested, material: PP, transparent, flat cap, suitable for qPCR, 12 piece(s)/minigrip bag |

### **Product characteristics**

|                        |  |
|------------------------|--|
| Cap                    | flat cap                                       |
| Closure type           | push cap                                       |
| Suitable for           | PCR plates (REF 72.1979.102/.132, 72.1982.202) |
| Suitable for qPCR      | yes  |
| Number of wells        | 8  |
| Autoclavable (max. °C) | 121 °C   |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

|                   |         |
|-------------------|---------|
| Width of product  | 7.4 mm  |
| Length of product | 80.5 mm |

### **Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | transparent        |

### **Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Purity standard  | PCR Performance Tested  |
| Batched          | yes                     |

### **Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1200          |
| Type of smallest subpackaging | minigrip bag  |
| Piece(s) / inner box          | 300           |
| Piece (s) / outer case        | 1200          |
| Piece(s) / pallet             | 151200        |
| Depth of case                 | 355 mm        |
| Width of case                 | 255 mm        |
| Height of case                | 123 mm        |
| Case volume                   | 0.0111 cbm    |
| Weight of product             | 0.0004 kg     |
| Weight of case                | 0.69 kg       |
| EAN of inner box              | 4038917336369 |
| EAN case                      | 4038917336352 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

