



### Product description

|                     |   |
|---------------------|---|
| Order number        | 72.1981   |
| Product description | PCR plate half skirt, 96 well, transparent, Low Profile, 100 µl, PCR Performance Tested, material: PP, 10 piece(s)/minigrip bag |

### Product characteristics

|                                 |               |
|---------------------------------|---------------|
| Suitable for qPCR               | yes           |
| Base shape                      | conical       |
| Construction of the plate edges | half skirt    |
| High or low profile             | Low Profile   |
| Number of wells                 | 96            |
| Volume of work                  | 5 µl - 100 µl |
| Centrifugation max (RCF)        | 3000 x g      |
| Centrifugation temperature      | 22            |
| Centrifugation time             | 5 Min         |
| Autoclavable (max. °C)          | 121 °C        |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

|                   |          |
|-------------------|----------|
| Sample volume     | 0.1 ml   |
| Max. volume       | 100 µl   |
| Diameter          | 5.48 mm  |
| Width of product  | 85.4 mm  |
| Height of product | 17.37 mm |
| Length of product | 127.8 mm |

### **Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | transparent        |

### **Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Purity standard  | PCR Performance Tested  |
| Batched          | yes                     |

### **Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 100           |
| Type of smallest subpackaging | minigrip bag  |
| Piece(s) / inner box          | 50            |
| Piece (s) / outer case        | 100           |
| Piece(s) / pallet             | 9800          |
| Depth of box                  | 270 mm        |
| Width of box                  | 230 mm        |
| Height of box                 | 96 mm         |
| Depth of case                 | 282 mm        |
| Width of case                 | 210 mm        |
| Height of case                | 247 mm        |
| Case volume                   | 0.0146 cbm    |
| Weight of product             | 0.014 kg      |
| Weight of case                | 1.75 kg       |
| EAN of inner box              | 4038917573474 |
| EAN case                      | 4038917573467 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

