



### **Product description**

|                     |  |
|---------------------|--|
| Order number        | 72.735.992   |
| Product description | PCR single tube, 0.5 ml, PCR Performance Tested, red/green/yellow/blue/violet/orange, material: PP, with anti-contamination shield, flat cap, per colour 500 piece(s)/minigrip bag |

### **Product characteristics**

|                            |                |
|----------------------------|----------------|
| Graduation                 | yes            |
| Cap                        | flat cap       |
| Closure type               | push cap       |
| Suitable for qPCR          | no             |
| Base shape                 | conical        |
| High or low profile        | High Profile   |
| Number of wells            | 1              |
| Volume of work             | 25 µl - 0,5 ml |
| Centrifugation max (RCF)   | 8000 x g       |
| Centrifugation temperature | 22             |
| Centrifugation time        | 5 Min          |
| Autoclavable (max. °C)     | 121 °C         |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**PCR single tube, 0.5 ml, PCR Performance**  
**Tested, colour mix, PP, flat cap**

**Page 2**

**Size**

|                   |         |
|-------------------|---------|
| Sample volume     | 0.5 ml  |
| Max. volume       | 500 µl  |
| Diameter          | 7.9 mm  |
| Width of product  | 10 mm   |
| Height of product | 29.5 mm |
| Length of product | 30.9 mm |

**Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | colour mix         |

**Purity & certification**

|                  |                                |
|------------------|--------------------------------|
| Product category | In vitro diagnostic device, CE |
| CE certified     | CE                             |
| Purity standard  | PCR Performance Tested         |
| Batched          | yes                            |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 3000          |
| Type of smallest subpackaging | minigrip bag  |
| Piece(s) / inner box          | 3000          |
| Piece (s) / outer case        | 3000          |
| Piece(s) / pallet             | 198000        |
| Depth of case                 | 428 mm        |
| Width of case                 | 323 mm        |
| Height of case                | 155 mm        |
| Case volume                   | 0.0214 cbm    |
| Weight of product             | 0.0004 kg     |
| Weight of case                | 1.67 kg       |
| EAN of inner box              | 4038917120067 |
| EAN case                      | 4038917017053 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-06-30

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**